



OSA Technology Sdn Bhd



1 MINI FRAGMENT

Regain Life
Live Life





OSA Technology Sdn Bhd

1 VISION STATEMENT

Regain Life .
Live Life

2 MISSION STATEMENT



2 To bring innovative ideas from surgeons to the market through systematic process flow from design, prototype, test, intellectual property and innovative marketing skill.

To introduce world class orthopaedics and medical devices to the global market.

To seek continuous improvement in all business areas and contribute to society via knowledge, skills and quality products and services.

3 ABOUT OSA TECHNOLOGY

3 OSA Technology Sdn Bhd (OSA) is an innovation led company providing quality and advance designed products to the dynamic global medical needs.



4 OSA conceptualises, designs and produces technologically advance medical products with cooperation from scientists in various fields utilizing the latest leading edge technology.

As an innovator and manufacturer, OSA quality standards are ensured by its Quality Management System. OSA embraces the relevant requirements of the International Organisation of Standard (I.S.O.) Code of Practices and EC Directive 93/42/EEC Annex II, Article 3 (CE).

4 OUR SPECIALITIES

OSA Regain Life . Live Life.

1 MINI FRAGMENT



Quarter Tubular Plate
2.7mm



Straight Mini Plate
2.0mm

APPL Code	Holes	Length
46.3020.01A	3	23mm
46.3020.01B	4	31mm
46.3020.01C	5	39mm
46.3020.01D	6	47mm
46.3020.01E	7	55mm
46.3020.01F	8	63mm

APPL Code	Holes	Length
46.3026.01A	3	17mm
46.3026.01B	4	23mm
46.3026.01C	5	29mm
46.3026.01D	6	35mm

REMARK | Implants not listed are available upon request with minimum order quantity.

1 MINI FRAGMENT



Self
Compression
Plates
2.7mm



Self
Compression
Plates
2.0mm

APPL Code	Holes	Length
46.3021.02A	2	20mm
46.3021.02B	2	26mm
46.3021.02C	3	28mm
46.3021.02D	4	36mm
46.3021.02E	5	44mm
46.3021.02F	6	52mm
46.3021.02G	7	60mm
46.3021.02H	8	68mm
46.3021.02I	9	76mm
46.3021.02J	10	84mm
46.3021.02K	11	92mm
46.3021.02L	12	100mm

APPL Code	Holes	Length
46.3021.01A	4	22mm
46.3021.01B	5	27mm
46.3021.01C	6	32mm
46.3021.01D	7	37mm
46.3021.01E	8	42mm

REMARK | Implants not listed is available upon request with minimum order quantity.

1 MINI FRAGMENT



Epiphysis Mini Plates
2.7mm




Epiphysis Mini Plates
2.0mm

APPL Code	Description
46.3011.02A	L-Oblique (Left)
46.3011.02B	L-Oblique (Right)
46.3011.02C	Left (90 Degree)
46.3011.02D	Right (90 Degree)
46.3011.02E	T - Plate

APPL Code	Description
46.3011.01A	L-Oblique (Left)
46.3011.01B	L-Oblique (Right)
46.3011.01C	Left (90 Degree)
46.3011.01D	Right (90 Degree)
46.3011.01E	T - Plate

REMARK | Implants not listed are available upon request with minimum order quantity.

1 MINI FRAGMENT



ø1.5mm
Cortical Screw
0.5mm Pitch
Full
Thread



ø2.0mm
Cortical Screw
0.6mm Pitch
Full
Thread

APPL Code	Length
46.0402.01A	6mm
46.0402.01B	7mm
46.0402.01C	8mm
46.0402.01D	9mm
46.0402.01E	10mm
46.0402.01F	11mm
46.0402.01G	12mm
46.0402.01H	14mm
46.0402.01I	16mm
46.0402.01J	18mm
46.0402.01K	20mm
46.0402.01L	22mm
46.0402.01M	24mm
46.0402.01N	26mm
46.0402.01O	28mm
46.0402.01P	30mm

APPL Code	Length
46.0402.02A	6mm
46.0402.02B	8mm
46.0402.02C	10mm
46.0402.02D	12mm
46.0402.02E	14mm
46.0402.02F	16mm
46.0402.02G	18mm
46.0402.02H	20mm
46.0402.02I	22mm
46.0402.02J	24mm
46.0402.02K	26mm
46.0402.02L	28mm
46.0402.02M	30mm
46.0402.02N	32mm
46.0402.02O	34mm
46.0402.02P	36mm
46.0402.02Q	38mm
46.0402.02R	40mm

REMARK | Implants not listed are available upon request with minimum order quantity.

1 MINI FRAGMENT



ø2.7mm
Cortical Screw
1.0mm Pitch
Full
Thread



ø2.7mm
Cortical Screw
1.0mm Pitch
Self-Tapping
Full
Thread

APPL Code	Length
46.0402.03A	8mm
46.0402.03B	10mm
46.0402.03C	12mm
46.0402.03D	14mm
46.0402.03E	16mm
46.0402.03F	18mm
46.0402.03G	20mm
46.0402.03H	22mm
46.0402.03I	24mm
46.0402.03J	26mm
46.0402.03K	28mm
46.0402.03L	30mm
46.0402.03M	32mm
46.0402.03N	34mm
46.0402.03O	36mm
46.0402.03P	38mm
46.0402.03Q	40mm
46.0402.03R	42mm
46.0402.03S	44mm
46.0402.03T	46mm
46.0402.03U	48mm
46.0402.03V	50mm

APPL Code	Length
46.0402.04A	6mm
46.0402.04B	8mm
46.0402.04C	10mm
46.0402.04D	12mm
46.0402.04E	14mm
46.0402.04F	16mm
46.0402.04G	18mm
46.0402.04H	20mm
46.0402.04I	22mm
46.0402.04J	24mm
46.0402.04K	26mm
46.0402.04L	28mm
46.0402.04M	30mm
46.0402.04N	32mm
46.0402.04O	34mm
46.0402.04P	36mm
46.0402.04Q	38mm
46.0402.04R	40mm

REMARK | Implants not listed are available upon request with minimum order quantity.



OSA Technology Sdn Bhd

Regain Life . **Live Life.**

A large, white, rectangular area with horizontal lines, resembling a sheet of lined paper. It is centered on the page and has a subtle drop shadow effect. A large, faint, grey watermark of the 'OSA' logo is visible in the center of this area.



Regain Life . **Live Life.**



OSA Technology Sdn Bhd

Regain Life . **Live Life.**

www.osatechnology.com
customercare@osatechnology.com

Corporate Office

OSA Technology Sdn.Bhd (643478-M)
OSA Niaga Sdn. Bhd (773898-V)
Lot 338, 3rd Floor, Block 5
Laman Seri Business Park
Persiaran Sukan, Seksyen 13
40100 Shah Alam
Selangor, Malaysia.

Tel : +603 5518 2335
Fax : +603 5510 2257

Manufacturing Plant

OSA Technology Sdn.Bhd (643478-M)
Block C, UKM-MTDC
Smart Technology Centre II
Universiti Kebangsaan Malaysia
43600 Bangi,
Selangor, Malaysia.

Tel : +603 8926 1351
Fax : +603 8926 1361